YUCERA

Dental Zirconia Block User Manual

| | | |
|---------------------------------|-----|-------|
| INSTRUCTIONS FOR USE | EN | 1-4 |
| INSTRUCTIONS POUR L'UTILISATION | FR | 5-8 |
| INSTRUCCIONES DE USO | ES | 9-12 |
| ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ | RU | 13-16 |
| INSTRUÇÕES DE USO | PT | 17-20 |
| مادختسالا تاميلعت | AR | 21-24 |
| ISTRUZIONI PER L'USO | IT | 25-28 |
| KULLANIM YÖNERGESİ | TUR | 29-32 |
| GEBRAUCHSANWEISUNG | DE | 33-36 |

CE 0197

EN Instructions for use

Manufacture Date

See label

Shenzhen Yurucheng Dental Materials Co., Ltd.

Manufacture Address: 101,201,301, Building A, No. 35, Zhuqing Rd., Shijing Community, Shijing Street, Pingshan District, Shenzhen.

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Tel: +31(0) 642758955 E-mail: ar@umedwings.eu

Please read the product manual before use, and make it by professional technicians; install, adjust and wear it by professional doctors; use it according to the instructions and doctor's instructions.

Date of formulation: 2023.07.11 Edition: A.4

Product Name Dental Zirconia Block

Organization Structure

Model White:

The main ingredients include zirconia, yttrium oxide other oxides.

$ZrO_2+HfO_2+Y_2O_3 \geq 99\%$

$Y_2O_3: 4.5\% \sim 6.0\%$

$Al_2O_3 \leq 0.5\%$

Other Oxides $\leq 0.5\%$

Model Color and Multilayer:

The main ingredients include zirconia, ferric oxide, Erbium oxide and other oxides.

$ZrO_2+HfO_2+Y_2O_3 \geq 98.2\%$

Erbium oxide $\leq 1.0\%$

Ferric oxide $\leq 0.3\%$

Other Oxides $\leq 0.5\%$

Service Life 5 years

Application method

Open the bag take out porcelain pieces, placed in special dental processing equipment for processing, after dyeing, high temperature (1450~ 1600 °C).

After sintering and crystallization, the finished product will be made after dressing, overlaying decorative porcelain or glazing and other processes. Finally, it will be used by professional doctors for the repair, installation and wearing of human denture.

Intended Use

The Dental Zirconia Block is used for making crowns, Bridges, inlays and veneers of fixed denture.

Physicochemical property

Density (post sintering): $\geq 6.0g/cm^3$

Class:5: Flexural strength (post sintering) $\geq 800MPa$

Class:4: Flexural strength (post sintering) $\geq 600MPa$

Chemical solubility: $< 100\mu g.cm^{-2}$

Radioactivity (uranium-238) $\leq 1.0Bq/g$

coefficient of linear thermal expansion: $(10.5 \pm 0.5) \times 10^{-6}K^{-1}$

Type: II, supplied in block

Class:5: White ST/White SHT/White HT/Color ST/Color SHT/Multilayer SHT/Multilayer ST.

Class:4: White UT/Color UT/Multilayer 3D/Multilayer 4D/Multilayer UT.

Notice

1. Please check well the package and inside product, if any damage, forbid to use.
2. Don't press and collide during store, transit and processing.
3. Unburned porcelain should not be directly used in human denture or oral repair.
4. This product needs to be made by professional technicians and used by professional doctors.
5. Expired product is forbidden to reuse.
6. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
7. Single use products, do not reuse.

Package Carton packaging.

Disposal Abandon directly, the products do not pollute the environment after waste.

Biological Properties

1. No cytotoxicity
2. No sensitization
3. No acute oral stimulation
4. No acute systemic toxicity
5. No chronic systemic toxicity
6. AMES test: negative

Storage

Storage at indoor temperature (generally 5 -50 °C), ventilated and dry conditions.

Contraindication

In the event that there is an insufficient occlusal clearance and/or vertical prep wall, making the preparation unsuitable for an all-ceramic restoration, an alternative material must be chosen. Inlay bridges, endosseous implants and root posts are other contraindications.

Suggested sintering curve

Sintering can take place in all common dental sintering furnaces which are approved for the sintering of restorations made of zirconium dioxide. Since zirconium dioxide is known to be a poor thermal conductor, it is recommended that you slowly heat the mounts to the required temperature (see sintering graph) and cool them down equally slowly.

Sintering Curve

HT/ST /ST color /ST Multilayer
Sintering Curve (1-5)

| Sintering step | Start Temperature(°C) | End Temperature(°C) | Time(Min) | Rate(°C/Min) |
|----------------|-----------------------|---------------------|-----------|--------------|
| step 1 | 20 | 1000 | 130 | 7.5 |
| step 2 | 1000 | 1530 | 156 | 3.4 |
| step 3 | 1530 | 1530 | 120 | 0 |
| step 4 | 1530 | 800 | 104 | -7 |
| step 5 | 800 | natural cooling 20 | / | / |

HT/ST /ST color /ST Multilayer
Sintering Curve (6-10)

| Sintering step | Start Temperature(°C) | End Temperature(°C) | Time(Min) | Rate(°C/Min) |
|----------------|-----------------------|---------------------|-----------|--------------|
| step 1 | 20 | 1000 | 245 | 4 |
| step 2 | 1000 | 1530 | 177 | 3 |
| step 3 | 1530 | 1530 | 120 | 0 |
| step 4 | 1530 | 800 | 146 | -5 |
| step 5 | 800 | natural cooling 20 | / | / |

SHT Multilayer/3D Plus Multilayer/UT Multilayer
4D Pro Multilayer Sintering Curve (1-5)

| Sintering step | Start Temperature(°C) | End Temperature(°C) | Time(Min) | Rate(°C/Min) |
|----------------|-----------------------|---------------------|-----------|--------------|
| step 1 | 20 | 1000 | 130 | 7.5 |
| step 2 | 1000 | 1500 | 147 | 3.4 |
| step 3 | 1500 | 1500 | 120 | 0 |
| step 4 | 1500 | 800 | 100 | -7 |
| step 5 | 800 | natural cooling 20 | / | / |

SHT Multilayer/3D Plus Multilayer/UT Multilayer
4D Pro Multilayer Sintering Curve (6-10)

| Sintering step | Start Temperature(°C) | End Temperature(°C) | Time(Min) | Rate(°C/Min) |
|----------------|-----------------------|---------------------|-----------|--------------|
| step 1 | 20 | 1000 | 245 | 4 |
| step 2 | 1000 | 1500 | 167 | 3 |
| step 3 | 1500 | 1500 | 120 | 0 |
| step 4 | 1500 | 800 | 140 | -5 |
| step 5 | 800 | natural cooling 20 | / | / |

HT/ST /ST color /ST Multilayer
Sintering Curve (11-14)

| Sintering step | Start Temperature(°C) | End Temperature(°C) | Time(Min) | Rate(°C/Min) |
|----------------|-----------------------|---------------------|-----------|--------------|
| step 1 | 20 | 1000 | 326 | 3 |
| step 2 | 1000 | 1530 | 265 | 2 |
| step 3 | 1530 | 1530 | 120 | 0 |
| step 4 | 1530 | 800 | 183 | -4 |
| step 5 | 800 | natural cooling 20 | / | / |

UT Sintering Curve (1-5)

| Sintering step | Start Temperature(°C) | End Temperature(°C) | Time(Min) | Rate(°C/Min) |
|----------------|-----------------------|---------------------|-----------|--------------|
| step 1 | 20 | 1000 | 130 | 7.5 |
| step 2 | 1000 | 1450 | 132 | 3.4 |
| step 3 | 1450 | 1450 | 120 | 0 |
| step 4 | 1450 | 800 | 93 | -7 |
| step 5 | 800 | natural cooling 20 | / | / |

SHT Multilayer/4D Pro Multilayer
Sintering Curve (11-14)

| Sintering step | Start Temperature(°C) | End Temperature(°C) | Time(Min) | Rate(°C/Min) |
|----------------|-----------------------|---------------------|-----------|--------------|
| step 1 | 20 | 1000 | 326 | 3 |
| step 2 | 1000 | 1500 | 250 | 2 |
| step 3 | 1500 | 1500 | 120 | 0 |
| step 4 | 1500 | 800 | 175 | -4 |
| step 5 | 800 | natural cooling 20 | / | / |

Wall and edge thickness

| Type | Anterior | Connector cross section | Posterior | Connector cross section |
|----------------|-----------|-------------------------|-----------|-------------------------|
| | Thickness | | Thickness | |
| Single crown | ≥0.9mm | ≥9mm ² | ≥1.0mm | ≥12mm ² |
| 3 unit bridges | ≥0.9mm | | | |
| Long bridges | ≥0.9mm | | | |

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